www.lagis.com.tw/eifu_en.php

BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

◆ DEVICE DESCRIPTION :

The LAGIS Insufflation Needle incorporates a spring-loaded blunt stylet mechanism for establishing pneumoperitoneum prior to trocar and cannula insertion in laparoscopic procedures. The Insufflation Needle consists of a 14-gauge stainless steel needle and spring loaded stylet that retracts backward as the needle is pushed into the abdominal tissue. When the needle tip punctures through the anterior wall of the abdomen, the stylet advances to expose gas port. For insufflation, the proximal end of the device is attached to a plastic handle which offers optimum control when inserting or removing from the abdomen during laparoscopic procedures. The LAGIS Insufflation Needle is equipped with a stopcock and a luer lock connector for attaching to an insufflation tubing.

Technical data:

Item	Tensile force	Flow rate	Insufflation pressure
Value	>30 kgf	1~3 L/min	5~20 mmHg

◆ INDICATIONS/INTENDED PURPOSE :

- The LAGIS Insufflation Needle is intended for providing an access for insufflating gases into body cavity during gynecological laparoscopy, laparoscopic cholecystectomy and other laparoscopic procedures.
- Expected clinical benefits: simple operation, little time-consumption and effectiveness, no incidence of major complications (life threatening) visceral and vascular injuries as compared to the open technology and direct trocar entry technique.
- Patient target group: LAGIS Insufflation Needle is intended for adults who need laparoscopic procedures.

◆ CONTRAINDICATIONS :

This device is not intended for use except as indicated. It is not intended for use when endoscopic techniques are contraindicated.

◆ INSTRUCTIONS FOR USE :

- Prior to use, inspect the contents of the packaging and the packaging itself to ensure the integrity is not compromised.
- Prior to insertion, inspect the tip of LAGIS Insufflation Needle to ensure the stylet retracts smoothly when a gentle pressure is applied to it. Upon release of the pressure, the stylet extends freely beyond the needle tip to its resting position. Do not attempt to use the device if the stylet does not function properly.
 - Perform peritoneal insertion using standard insertion technique.
 - Advance the needle to desired insertion depth.
 - Attach the proximal end to a tubing connector.
 - Start insufflating.
 - Withdraw the Insufflation Needle when desired pressure is achieved or when the laparoscopic procedure is completed.
 - Dispose of the Insufflation Needle properly and safely by healthcare professionals in accordance with the national waste management regulation.

◆ WARNINGS AND PRECAUTIONS :

- Do not use if the device or the sterile packaging is open or damaged, as this may cause an increase in the incidence of wound infections.
- If endoscopic instruments or accessories from different manufacturers are used together in a procedure, verify compatibility prior to performing the procedures.
- This device should be operated only by surgeons with adequate training and familiarity with surgical techniques. Consult medical literature for further information regarding techniques, hazards, contra-indications, and complications prior to performing the procedures.
- This device is packaged and sterilized for single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocess, or reesterilization may create contamination, infection, cross-infection, or even compromise the integrity of the device and lead to device failure which in turn may result in patient injury, illness or death.

- Do not attempt to use the device if the stylet is not functioning properly, such as safety indicator does not move back together, it means needle tip cannot be exposed for abdominal wall puncture.
- Return the device and packaging to the distributor or implement safe disposal handling by the healthcare professionals if the sterile packaging is damaged or unintentionally opened before use.
- The used device which comes into contact with bodily fluids should be considered as biomedical waste and safely disposed of by healthcare professionals in compliance with the national waste management regulation to prevent patients, users, and other persons from biological contamination or infection.
- In the event of malfunction of the device, changes in performance that they may affect safety, or even any serious incident, please report to the manufacturer, the competent authority of the Member State and the authorized representative immediately according to the contact information provided on the label or in this instructions for use.
- This product and its packaging have been sterilized with ethylene oxide. The packaging may expose the user and patient to ethylene oxide, a chemical known to cause cancer or birth defects or other reproductive harm.
- Improper assembly or operation of the device may result in prolong surgical time, and tissue or organ injuries from mechanical damage to the patient or operator.

Compatibility:

Insufflation tubing connectors and CO₂ Insufflator.

◆ ENVIRONMENTAL CONDITIONS FOR STORAGE :

Appropriate storage environment is a clean and dry area away from sunlight with a temperature range of 13~30°C (55.4~86°F).

◆ EXPLANATION OF SYMBOLS :

	Temperature limit		Caution		Authorized representative in the European Community/European Union
	Do not re-sterilize		Manufacturer		Not made with natural rubber latex
	Do not re-use		Consult instructions for use or consult electronic instructions for use		Do not use if package is damaged and consult instructions for use
	Date of manufacture		Single sterile barrier system		Sterilized using ethylene oxide
	Medical device		Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.		

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