www.lagis.com.tw/eifu_en.php

BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

◆ DEVICE DESCRIPTION :

The LAGIS Irrigation Pump is intended to pressurize an irrigation solution bag to support high pressure or high flow of irrigation solutions during surgery. It utilizes pressurized gas to provide irrigation pressure up to 700 mmHg. The user controls the flow rate using a roller clamp on the irrigation set. The maximum flow rate is determined by the viscosity of the fluid, the characteristics of the irrigation set, and the pressure applied.

PRODUCT MODEL

IP-1000

IP-2000

◆ INDICATIONS/INTENDED PURPOSE :


- The LAGIS Irrigation Pump has application in general surgical procedures for rapid delivery of irrigation solutions where flow rate precision is not required.
- Expected clinical benefits: The device applied to clear debris or blood can be helpful to clearly identify the origin site of bleeding during an operative procedure. It may also reduce the rate of surgical site infections through the removal of dead or damaged tissue, metabolic waste, and wound exudate.
- Patient target group: the LAGIS Irrigation Pump is intended for adults under the surgery that requires a rapid administration of fluid products.

◆ CONTRAINDICATIONS :

The device is not intended for infusion of blood, IVs, cystoscopy, or for hysteroscopy.

◆ ILLUSTRATION OF IRRIGATION PUMP COMPONENTS

1. Control Switch
2. Pressure Regulator Knob
3. Rapid Deflation Switch
4. Latch
5. Cover
6. Air Source Connector
7. Air Bladder (behind Cover)
8. Bag Hanger
9. Spring
10. Pressure Gauge

◆ INSTRUCTIONS FOR USE :

SPECIFICATION OF THE LAGIS IRRIGATION PUMP

Input	Compressed air, nitrogen (Do not exceed 60 psi.)
Gauge accuracy	± 30 mmHg
Operating Pressure	0-700 mmHg
Compatible Fluid Bags	Flexible bags, supplied by manufacturers are compatible with Irrigation Pump. Rigid and semi-rigid containers are not compatible.

OPERATING INSTRUCTIONS

1. Fix Irrigation Pump to IV pole outside the sterile field.
2. Flip Control Switch to "DEPRESSURIZE" before opening the Cover.
3. Connect the Air Source Connector to the compressed nitrogen or air bottles/tanks for applying pressure to the Irrigation Pump.
Note: Connecting to a source other than those recommended, including the fluid warming system, may result in incorrect operating pressures or damage to internal components. Do not connect to nitrous oxide or other anesthetic gases. There is no indication on the Irrigation Pump that a gas pressure source has been connected.
4. Open the Cover and hang the fluid bag on the Bag Hanger.
5. Spike the fluid bag with an irrigation set. Use the irrigation set per the instructions for use of the irrigation set.

6. Close the Cover.
Note: If the air bladder is partly inflated, turn on the rapid deflation switch to deflate the air bladder. It will exhaust residual air from the air bladder so the cover can be closed and latched.
7. Secure the Latch.
Note: Ensure the cover is latched before pressurizing.
8. Flip the Control Switch to "PRESSURIZE" and adjust the pressure by pulling the Pressure Regulator Knob and rotating clockwise to increase, and counter-clockwise to decrease. Push the Pressure Regulator Knob to lock the pressure.
Note: Do not operate the LAGIS Irrigation Pump at a pressure of 700 mmHg. The pressure has been set at the factory to approximately 300 mmHg.
9. Adjust the roller clamp on the irrigation set to achieve the desired flow rate.
Note: The LAGIS Irrigation Pump controls the pressure applied to the fluid bag. There are no alarms provided for empty bags, occlusions, or excessive flow. It is the user's responsibility to ensure proper fluid flow, including monitoring for excessive flow or occlusion.
10. Change the fluid bags: Flip Control Switch to "DEPRESSURIZE" and Rapid Deflation Switch to "ON". After the Air Bladder is deflated, open the Cover and replace the fluid bags.

CLEANING

The exterior of the pump should be wiped down with a damp cloth and germicide. Failure of the door and other parts may be caused by the application of cleaning solutions containing alcohol or solvents. No cleaning is required inside the pump. The use of other cleaning agents may result in premature failure of the product. Do not immerse the Irrigation Pump in liquids or use alcohol or solvents to clean. Do not autoclave.

PREVENTIVE MAINTENANCE

- Monthly, the user should:
 1. Examine the surface for cracks that may occur (e.g. due to the use of improper cleaning solutions). If cracks are observed, remove the Irrigation Pump from service and replace the damaged part(s).
 2. Check gauge for zero reading (within one division) with no pressure applied. A non-zero reading is an indication that Irrigation Pump may have been dropped or otherwise mishandled and the unit should be thoroughly checked.
 3. Check the regulator operating range of 0 to 700 mmHg. If the full pressure cannot be achieved, check for air leaks in tubing or air bladder.
- Annually, in addition to monthly checks, it is recommended to users:
 1. Verify maximum pressure is between 600 and 700 mmHg with 60 psi input.
 2. Check tubing and connections for leaks and signs of wear.
Note: The calibration cannot be performed by users. Factory calibration and replacement components are available if required.
 3. The user should replace the cover and air bladder if any cracks or damages are observed.

DISPOSAL

The primary plastic materials used in the Irrigation Pump are Polycarbonate and Nylon. The method of disposal of the Irrigation Pump at the end of life should be appropriate for these materials according to the relevant national laws.

◆ WARNINGS AND PRECAUTIONS :

- The storage period shown on the label is the period of the entire equipment. The air bladder should maintain its functionality for use 10000 times, and the air bladder will not deflate when the times limit is reached. If there is a problem with the air bladder or other components, please refer to the TROUBLESHOOTING section or contact the manufacturer for the service required.
- This device should be operated only by surgeons with adequate training and familiarity with surgical techniques. Consult medical literature for further information regarding techniques, hazards, contra-indications, and complications prior to performing the procedures.
- Improper assembly or operation of the device may result in prolonged surgical time and undesirable tissue injuries to the patient.
- To avoid the patient or user hurt, do not use the device when it is damaged or there is any component loosening.
- The use of excess irrigation fluid and irrigation time can lead to electrolyte imbalance, which is one of the complications, as this may cause prolonged surgical time and unintended hurt to the patient.
- High-pressure irrigation also used for hydro-dissection under an operative procedure may result in the patient hurt.
- If air sources or accessories from different manufacturers are used together with the LAGIS Irrigation Pump, verify compatibility prior to performing the procedure.
- In the event of malfunction of the device, changes in performance that may affect safety, or even any serious incident, please report to the manufacturer, the competent authority of the Member State and the authorized representative immediately according to the contact information provided on the label or in the instructions for use.
- The device are provided non-sterile, reusable and are to be used outside the sterile field only. The device must be thoroughly cleaned before the first use and after every use.


TROUBLESHOOTING

Problem	Cause	Solution
	Pump not hooked up to gas source	Connect pump to gas source
No Irrigation	Unit is not "ON"	Flip control switch to "PRESSURIZE"
	Irrigation pressure is set to zero	Pull knob and turn clockwise to reach desired pressure. Push to lock
	Rapid deflation switch is "ON"	Flip rapid deflation switch to "OFF"
	Kinked tubing	Straighten or replace tubing
	Fluid bag is empty	Replace fluid bag
	Fluid line clamp(s) is closed	Open clamp(s)
	Air bladder or internal air leak	Service required
Air Bladder not deflate	Control switch is not "OFF"	Flip control switch to "DEPRESSURIZE"
	Rapid deflation Switch is "OFF"	Flip rapid deflation switch to "ON"
	Switch damaged	Service required
	Internal problem	Service required
Leakage of solution	Bag not spiked securely	Secure spike in fluid bag
Unable to adjust the pressure	Incorrect pressure at the air source	Adjust the pressure of air source to 60 psi
Unable to reduce pressure	Damaged regulation knob or pressure gauge	Service required
	Unit is not "OFF"	Turn control switch to "DEPRESSURIZE"

◆ ENVIRONMENTAL CONDITIONS FOR STORAGE :

Appropriate storage environment is a clean and dry area away from sunlight with a temperature range of 13~30°C (55.4~86°F).

◆ EXPLANATION OF SYMBOLS :

	Not made with natural rubber latex		Manufacturer		Consult instructions for use (or consult electronic instructions for use).
	Caution		Serial Number		Federal(USA) law restricts this device to sale by or on the order of a physician.
	Temperature limit		Medical device		Authorized representative in the European Community / European Union

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**LAGIS[®] ENTERPRISE CO., LTD.**