www.lagis.com.tw/eifu\_en.php

**BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.**

## ◆ DEVICE DESCRIPTION:

The LAGIS Endo Bag is a single use device sterilized by ethylene oxide. It is composed of a string, an introducer tube, and an elastic plastic bag. The introducer tube allows the elastic plastic bag to deploy tissue/specimen collection in the body cavity through a proper size trocar. The string facilitates closure of the bag once the tissue/specimen has been collected.

## PRODUCT MODEL

The final product codes for various size Endo Bags are listed as follows:

Device Name	Product Code	Descriptions	Total volume (ml)	Pouch diameter (mm)
Endo Bag	EB-0507	5" x 7" specimen retrieval bag	300	80
Endo Bag	EB-0709	7" x 9" specimen retrieval bag	1000	110

## ◆ INDICATIONS/INTENDED PURPOSE:

- The LAGIS Endo Bag is intended for safe removal or collection of tissue specimens from the abdominal cavity without contamination and infection during laparoscopic surgical procedures.
- Expected clinical benefits: avoid accidental spillage of content or occult malignant cells in the abdominal cavity and contamination to abdominal cavity and wound tract.
- Patient target group: the LAGIS Endo Bag is intended for adults who need laparoscopic surgical procedures.

## ◆ CONTRAINDICATIONS:

- The LAGIS Endo Bag is not intended for use when minimally invasive techniques are contraindicated.
- The LAGIS Endo Bag is not intended for use with any tissue that may not fit within the confines of the bag and allow complete closure of the bag.

## ◆ INSTRUCTIONS FOR USE:

Prior to use, inspect the contents of the packaging and the packaging itself to ensure the integrity is not compromised.

- Leave the string accessible outside the body cavity and insert the LAGIS Endo Bag through a 10 mm or larger size trocar.
- Use a blunt tip endoscopic instrument to grasp the bag out of the introducer tube and unroll the bag gently to prevent any damage. Withdraw the introducer tube from the trocar after the bag begins to deploy.
- Place the specimen/tissue into the bag. Ensure the entire specimen/tissue fits within the bag.
- Once the specimen/tissue is completely within the bag, pull the string until the bag is tightly closed and leave it in the body cavity for removal at the end of the procedure.
- Remove the trocar beforehand and pull up the string to withdraw the bag through the incision. If the bag with the specimen/tissue cannot be withdrawn, carefully enlarge the incision for easier removal.

## ◆ WARNINGS AND PRECAUTIONS:

- Do not use if the device or the sterile packaging is open or damaged, as this may cause an increase in the incidence of wound infections.
- Do not pull up the bag through the trocar sleeve while the specimen/tissue is within the bag.
- If trocars or accessories from different manufacturers are used together in a procedure, verify compatibility prior to performing the procedure.
- This device is packaged and sterilized for single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocess, or resterilization may create contamination, infection, cross-infection, or even compromise the integrity of the device and lead to device failure which in turn may result in patient injury, illness or death.
- This device should be operated only by surgeons with adequate training and familiarity with surgical techniques. Consult medical literature for further information regarding techniques, hazards, contra-indications, and complications prior to performing the procedures.
- Return the device and packaging to the distributor or implement safe disposal handling by the healthcare professionals if the sterile packaging is damaged or unintentionally opened before use.
- The used device which comes into contact with bodily fluids should be considered as biomedical waste and safely disposed of by healthcare professionals in compliance with the national waste management regulation to prevent patients, users, and other persons from biological contamination or infection.
- In the event of malfunction of the device, changes in performance that they may affect safety, or even any serious incident, please report to the manufacturer, the competent authority of the Member State and the authorized representative immediately according to the contact information provided on the label or in this instructions for use.
- Improper assembly or operation of the device may result in prolong surgical time, patient allergy, and tissue or organ injuries from mechanical damage to the patient or operator.

## ◆ ENVIRONMENTAL CONDITIONS FOR STORAGE:

Appropriate storage environment is a clean and dry area away from sunlight with a temperature range of 13~30°C (55.4~86°F).

## ◆ EXPLANATION OF SYMBOLS:

	Temperature limit		Caution		Authorized representative in the European Community/European Union
	Do not resterilize		Manufacturer		Not made with natural rubber latex
	Do not re-use		Consult instructions for use or consult electronic instructions for use		Do not use if package is damaged and consult instructions for use
	Date of manufacture		Single sterile barrier system		Sterilized using ethylene oxide
	Medical device		Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.		

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