
www.lagis.com.tw/eifu\_en.php

**BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.**

## ◆ DEVICE DESCRIPTION :

The LAGIS Specimen Retrieval Bag is a single use device sterilized by gamma irradiation. It is composed of an elastic plastic bag, a push/pull plunger, a string knotting around the ring, and an introducer tube. The push/pull plunger and the introducer tube facilitates bag opening and tissue/specimen collection in the body cavity through a proper size trocar. In the fully deployed condition, the bag opening is maintained in a fully-open position by a metallic rim. The string enables the bag to close by pulling the ring after the tissue/specimen collection is completed.

- A: Metallic Rim
- B: Specimen Bag
- C: Introducer Tube
- D: Handle
- E: Blade Slot
- F: Plunger
- G: String
- H: Ring

## ◆ INDICATIONS/INTENDED PURPOSE :

- The LAGIS Specimen Retrieval Bag is intended for safe removal or collection of tissue specimens from the abdominal cavity without contamination and infection during laparoscopic surgical procedures.
- Expected clinical benefits: avoid accidental spillage of content or occult malignant cells in the abdominal cavity and contamination to abdominal cavity and wound tract.
- Patient target group: the LAGIS Specimen Retrieval Bag is intended for adults who need laparoscopic surgical procedures.

## ◆ CONTRAINDICATIONS :

- The LAGIS Specimen Retrieval Bag is not intended for use when minimally invasive techniques are contraindicated.
- The LAGIS Specimen Retrieval Bag is not intended for use with any tissue that may not fit within the confines of the bag and allow complete closure of the bag.
- Do not use morcellators with the LAGIS Specimen Retrieval Bag.

## ◆ INSTRUCTIONS FOR USE :

Prior to use, inspect the contents of the packaging and the packaging itself to ensure the integrity is not compromised.

- Insert the LAGIS Specimen Retrieval Bag through a 10 mm or larger size trocar. Ensure the printing "THIS SIDE UP" is facing up, which indicates that the mouth of the bag will be in the proper position when deployed.
- Push the plunger forward until it touches the handle to facilitate the bag opening. The bag may be fully unrolled by using a blunt tip endoscopic instrument.
- Place the specimen/tissue into the bag. Ensure the entire specimen/tissue fits within the bag.
- Once the specimen/tissue is completely within the bag, rotate the ring knotted on the string away from the plunger.
- Pull the ring until the bag is tightly closed and detached from the metallic rim. Use a blunt tip endoscopic instrument to aid detaching if needed.
- Cut the string by the blade slot located on the handle.
- Pull the plunger away from the handle until the metallic rim is retracted completely into the introducer tube.
- Withdraw the LAGIS Specimen Retrieval Bag from the trocar and leave the bag with specimen/tissue in the body cavity for removal at the end of the procedure.

- Remove the trocar beforehand and pull up the string to withdraw the bag through the incision. If the bag with the specimen/tissue cannot be withdrawn, carefully enlarge the incision for easier removal.

## ◆ WARNINGS AND PRECAUTIONS :

- Do not use if the device or the sterile packaging is open or damaged, as this may cause an increase in the incidence of wound infections.
- Do not pull up the bag through the trocar sleeve while the specimen/tissue is within the bag.
- Do not push the plunger to facilitate the bag opening until the tip of the introducer tube is passed through the trocar sleeve.
- Prior to withdrawing the LAGIS Specimen Retrieval Bag through the trocar sleeve, ensure the metallic rim is fully retracted into the introducer tube.
- If trocars or accessories from different manufacturers are used together in a procedure, verify compatibility prior to performing the procedure.
- This device is packaged and sterilized for single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocess, or resterilization may create contamination, infection, cross-infection, or even compromise the integrity of the device and lead to device failure which in turn may result in patient injury, illness or death.
- This device should be operated only by surgeons with adequate training and familiarity with surgical techniques. Consult medical literature for further information regarding techniques, hazards, contra-indications, and complications prior to performing the procedures.
- Return the device and packaging to the distributor or implement safe disposal handling by the healthcare professionals if the sterile packaging is damaged or unintentionally opened before use.
- The used device which comes into contact with bodily fluids should be considered as biomedical waste and safely disposed of by healthcare professionals in compliance with the national waste management regulation to prevent patients, users, and other persons from biological contamination or infection.
- Do not attempt to repair the device. The use of a defective device or a device that becomes defective during surgery may result in loose parts that could drop into the field of operation.
- In the event of malfunction of the device, changes in performance that they may affect safety, or even any serious incident, please report to the manufacturer, the competent authority of the Member State and the authorized representative immediately according to the contact information provided on the label or in this instructions for use.
- Improper assembly or operation of the device may result in prolong surgical time, patient allergy, and tissue or organ injuries from mechanical damage to the patient or operator.
- If the endo bag containing the specimen cannot be removed through the puncture site, the puncture site must be carefully enlarged to facilitate the removal of the bag. Do not force the endo bag taken out without enlarging the puncture site. If so, this may cause the endo bag to rupture, and the contents spill out of the bag.

## ◆ ENVIRONMENTAL CONDITIONS FOR STORAGE :

Appropriate storage environment is a clean and dry area away from sunlight with a temperature range of 13~30°C (55.4~86°F).

## ◆ EXPLANATION OF SYMBOLS :

|  |                     |  |   |  |  |
|--|---------------------|--|---|--|--|
|  | Temperature limit   |  | Caution   |  | Authorized representative in the European Community/European Union |
|  | Do not resterilize  |  | Manufacturer  |  | Not made with natural rubber latex                                 |
|  | Do not re-use       |  | Consult instructions for use or consult electronic instructions for use                     |  | Do not use if package is damaged and consult instructions for use  |
|  | Date of manufacture |  | Single sterile barrier system   |  | Sterilized using irradiation                                       |
|  | Medical device      |  | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. |  |  |

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