www.lagis.com.tw/eifu\_en.php

**BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.**

**◆ DEVICE DESCRIPTION :**

The LAGIS Mini Laparoscopy Grasping Device is a single use device sterilized by ethylene oxide. It consists of a suture passer and cones in different diameter. The suture passer has a shaft with a notched tip which is used to capture and hold the suture and a syringe-like handle that facilitates the opening/retracting of the tip by depressing the button. The cones which come in 5 mm, 12 mm, and 15 mm have holes in the cylindrical head to allow the shaft to pass through it.


**1. "LAGIS" Mini Laparoscopy Grasping Device**


**2. 5mm cone**

**3. 12mm cone**

**4. 15mm cone**

**◆ INDICATIONS/INTENDED PURPOSE :**

- The LAGIS Mini Laparoscopy Grasping Device has applications in laparoscopic and endoscopic procedures to facilitate percutaneous suturing for closing incision sites and fixation of the surgical mesh.
- Expected clinical benefits: reduce the surgery risk, shorten the surgery time, minimize the wound, reduce the bleeding situation and easily controlled for the clinician.
- Patient target group: The LAGIS Mini Laparoscopy Grasping Device is intended for adults who need laparoscopic and endoscopic procedures.

**◆ CONTRAINDICATIONS :**


The LAGIS Mini Laparoscopy Grasping Device is not intended for use when it is not possible to visualize the device during insertion and removal or when laparoscopic techniques are contraindicated.

**◆ INSTRUCTIONS FOR USE :**

Prior to use, inspect the contents of the packaging and the packaging itself to ensure the integrity is not compromised.

**Closure of fascial trocar wound: (MG-10512N / MG-11215N: Suture passer with cones)**


1. Once the procedure is completed, withdraw the trocar and insert the appropriately sized cone into the incision. Under the laparoscopic visibility, ensure the cone is perpendicular to the peritoneum.


2. Open the notched tip by depressing the button and place the suture (size 0, 2-0, 3-0) into it.


3. Release the button to close the notched tip on the suture. With the tip in the closed position, insert the LAGIS Mini Laparoscopy Grasping Device through the cone into the fascia and peritoneum. Drop the suture and withdraw the LAGIS Mini Laparoscopy Grasping Device.


4. With the notched tip in the closed position, insert the LAGIS Mini Laparoscopy Grasping Device through the opposite of the cone to pick up the suture. Pull the suture up through the cone.


5. Once the suture is in place, remove the cone and tie the suture to complete the incision closure.  
Note: Decrease the abdominal pressure to reduce the distention for better approximation of the incision site.


6. If using the 15 mm cone, insert the cone with hole 1 and 2 (refer to the figure below) aligned perpendicular to the incision. Perform steps 2) to 4) by using the holes labeled 1 and 2.

7. With another piece of suture, repeat steps 2) to 4) by using holes 3 and 4.


8. Once the sutures are in place, remove the cone and tie each suture to complete the incision closure.


**Hernia repair: (MG-10512N / MG-11215N: Suture passer)**

1. Measure the size of the wound and trim the surgical mesh according to the surgeon's discretion.
2. Thread four non-absorbable sutures (size 0, 2-0, 3-0), through the angles of the surgical mesh equidistantly and tie a knot.

3. Roll the surgical mesh lengthwise and introduce it into the abdomen through the trocar.
4. Make four small incisions corresponding to the angles of the surgical mesh to allow the LAGIS Mini Laparoscopy Grasping Device to insert through the fascia, muscle and peritoneum and into the abdomen.

5. Pick up the suture and pull it up through the fascia, muscle and peritoneum.
6. Repeat steps 4) and 5) with the other three tied sutures.

7. Secure the surgical mesh with staples/tacks. It is suggested that the staples/tacks should be placed in a proper distance from the edge of the surgical mesh for the best fixation points.

### ◆ WARNINGS AND PRECAUTIONS :

- Do not use if the device or the sterile packaging is open or damaged, as this may cause an increase in the incidence of wound infections.
- The pointed tip of the LAGIS Mini Laparoscopy Grasping Device may injure the tissues of the patients or users and other persons in contact with the tip. The pointed tip should be covered by the plastic protector at all times when the device is not in use to prevent unintentional injuries.
- Do not depress the button during withdrawal of the device, or it may result in inadvertent grasp/damage of the tissue.
- This device is packaged and sterilized for single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocess, or reesterilization may create contamination, infection, cross-infection, or even compromise the integrity of the device and lead to device failure which in turn may result in patient injury, illness or death.
- This device should be operated only by surgeons with adequate training and familiarity with surgical techniques. Consult medical literature for further information regarding techniques, hazards, contra-indications, and complications prior to performing the procedures.
- Return the device and packaging to the distributor or implement safe disposal handling by the healthcare professionals if the sterile packaging is damaged or unintentionally opened before use.
- The used device which comes into contact with bodily fluids should be considered as biomedical waste and safely disposed of by healthcare professionals in compliance with the national waste management regulation to prevent patients, users, and other persons from biological contamination or infection.
- Do not attempt to repair the device. The use of a defective device or a device that becomes defective during surgery may result in loose parts that could drop into the field of operation.
- In the event of malfunction of the device, changes in performance that they may affect safety, or even any serious incident, please report to the manufacturer, the competent authority of the Member State and the authorized representative immediately according to the contact information provided on the label or in this instructions for use.
- Improper assembly or operation of the device may result in tissue or organ injuries from mechanical damage to the patient or operator.

### ◆ ENVIRONMENTAL CONDITIONS FOR STORAGE :

Appropriate storage environment is a clean and dry area away from sunlight with a temperature range of 13~30°C (55.4~86°F).

### ◆ EXPLANATION OF SYMBOLS :

|  |                     |  |   |  |  |
|--|---------------------|--|---|--|--|
|  | Temperature limit   |  | Caution   |  | Authorized representative in the European Community/European Union |
|  | Do not re-sterilize |  | Manufacturer  |  | Not made with natural rubber latex                                 |
|  | Do not re-use       |  | Consult instructions for use or consult electronic instructions for use                     |  | Do not use if package is damaged and consult instructions for use  |
|  | Date of manufacture |  | Single sterile barrier system   |  | Sterilized using ethylene oxide                                    |
|  | Medical device      |  | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. |  |  |

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