[www.lagis.com.tw/eifu\\_en.php](http://www.lagis.com.tw/eifu_en.php)

**BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.**
**◆ DEVICE DESCRIPTION :**

The LAGIS Endoscopic Instruments- Scissors are comprised of a handle, a rotation knob, an insulated shaft, a pair of tips, and a cautery connector. The rotation knob located on the handle rotates the shaft 360°. To reposition the tips, turn the rotation knob. The tips are closed and opened by squeezing and releasing the handles. The cautery connection located on the top of the handle may be used for unipolar cautery when attached to standard cautery cables and their generators.

**PRODUCT MODEL**

| Scissors' Style | Insulated Shaft |         |             |
|-----------------|-----------------|---------|-------------|
|                 | Code            | OD (mm) | Length (mm) |
| Curved scissors | EI-330S         | 5       | 330         |

**◆ INDICATIONS/INTENDED PURPOSE :**

- The LAGIS Endoscopic Instruments- Scissors have applications in minimally invasive procedures for cutting and cauterization of tissue.
- Expected clinical benefits: Reduce the surgery risk, shorten the surgery time, minimize the wound, reduce the bleeding situation, and easily controlled for the clinician.
- Patient target group: The LAGIS Endoscopic Instruments-Scissors are intended for adults who need minimally invasive procedures.

**◆ CONTRAINDICATIONS :**

DO NOT USE the LAGIS Endoscopic Instruments - Scissors for contraceptive coagulation of fallopian tissue or when minimally invasive techniques are contraindicated; however, they may be used to achieve hemostasis following transection of the fallopian tube.

**◆ INSTRUCTIONS FOR USE :**

- Prior to use, inspect the contents of the packaging and the packaging itself to ensure the integrity is not compromised. The insulating layer of the shaft shall be free from any signs of cracks. All connectors and handles shall not be loose or wrongly connected.
- When attached to a high-frequency generator, make sure that the connectors are securely attached, and the HF power output is within the specified range.
  - LAGIS Endoscopic Instruments-Scissors can be used only when the abdominal cavity is inflated, and the laparoscope is already set up. Insert the device through a cannula into the body cavity to perform endoscopic surgical operations.
  - Upon completion of the procedure, retract LAGIS Endoscopic Instruments-Scissors and dispose of them properly and safely by healthcare professionals in accordance with the national waste management regulation.

**◆ WARNINGS AND PRECAUTIONS :**

- Do not use if the device or the sterile packaging is open or damaged, as this may cause an increase in the incidence of wound infections.
- If trocars or accessories from different manufacturers are used together in the procedure, verify compatibility prior to performing the procedure and ensure that electrical insulation or grounding is not compromised.

- This device is packaged and sterilized for single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocess, or resterilization may create contamination, infection, cross-infection, or even compromise the integrity of the device and lead to device failure which in turn may result in patient injury, illness or death.
- This device should be operated only by surgeons with adequate training and familiarity with surgical techniques. Consult medical literature for further information regarding techniques, hazards, contra-indications, and complications prior to performing the procedures.
- Improper assembly or operation of the device may result in tissue or organ injuries from mechanical or thermal damage to the patient or operator.
- To avoid the possibility of foreign bodies accidentally left in the patient's body, do not use the device when it is damaged or there is any component loosening.
- If failing to achieve the desired electrosurgical effect with normal power level settings, check the grounding pad ("return electrode") is properly affixed to the correct position and that all electrosurgical cables are firmly connected to the electrosurgical unit (ESU).
- Return the device and packaging to the distributor or implement safe disposal handling by the healthcare professionals if the sterile packaging is damaged or unintentionally opened before use.
- The used device which comes into contact with bodily fluids should be considered as biomedical waste and safely disposed of by healthcare professionals in compliance with the national waste management regulation to prevent patients, users, and other persons from biological contamination or infection.
- Do not attempt to repair the device. The use of a defective device or a device that becomes defective during surgery may result in loose parts that could drop into the field of operation.
- Do not use power generators that may deliver more than 2600 Volts peak to the Endoscopic Instruments. The electrosurgical generator power level should be set as low as necessary to achieve the desired electrosurgical effect. Excessive power levels may result in instrument malfunction, insulation failure, and injury or thermal damage to patients or users.
- Always use the lowest possible output setting necessary to achieve the desired electrosurgical effect. In order to lessen the possibility of unintended burn injury, use the active electrode only for the minimum time necessary. Higher settings may result in spurious insulation failure alarms and/or insulation breakdown.
- In the event of malfunction of the device, changes in performance that they may affect safety, or even any serious incident, please report to the manufacturer, the competent authority of the Member State and the authorized representative immediately according to the contact information provided on the label or in this instructions for use.

Compatibility:

4mm pin monopolar connectors and high-frequency generator.

**◆ ENVIRONMENTAL CONDITIONS FOR STORAGE :**

Appropriate storage environment is a clean and dry area away from sunlight with a temperature range of 13~30°C (55.4~86°F).

**◆ EXPLANATION OF SYMBOLS :**

|  |                     |  |   |  |   |
|--|---------------------|--|---|--|---|
|  | Temperature limit   |  | Caution   |  | Authorized representative in the European Community/European Union                          |
|  | Do not resterilize  |  | Manufacturer  |  | Not made with natural rubber latex  |
|  | Do not re-use       |  | Consult instructions for use or consult electronic instructions for use |  | Do not use if package is damaged and consult instructions for use                           |
|  | Date of manufacture |  | Single sterile barrier system   |  | Sterilized using ethylene oxide   |
|  | Medical device      |  | Type BF applied part  |  | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. |

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