www.lagis.com.tw/eifu_en.php

BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

◆ **DEVICE DESCRIPTION:**

The LAGIS Trocar is a single-use device sterilized by ethylene oxide. It is composed of a cannula, an obturator and a stopcock. The diameter of the cannula ranging from 3mm-12mm accommodates different sizes of endoscopic instruments. In order to be used in different patients and surgery, the length of the device comes in three sizes: Short, Standard and Long. The surface of the cannula also comes in two types: with thread and without thread.

There are five types of obturators: shielded tip, conical tip, blunt tip, dilating tip and optical tip, which offers a variety of configurations to meet users' needs.

- LAGIS Surgical Shielded Tip Trocar has a bladed tip with a shield designed to cover the cutting blade when the tip enters the body cavity.

- LAGIS Surgical Conical Tip Trocar has a rigid conical tip, which is sharp enough to enter the body cavity.

- LAGIS Surgical Blunt Tip Trocar has a blunt tip designed to prevent additional penetration of internal tissue.

- LAGIS Surgical Dilating Tip Trocar has a bladeless dilating tip designed to minimize abdominal wall and vessel trauma.

- LAGIS Surgical Optical Tip Trocar has a clear and transparent tip designed to offer direct visualization with the endoscope during insertion.

The cannula contains a set of seals to minimize the risk of air leakage when endoscopic instruments are inserted or withdrawn through the LAGIS Trocar. A stopcock is located at the side of the cannula for the connection to the luer lock to provide gas insufflation and desufflation.

◆ **INDICATIONS/INTENDED PURPOSE:**

- The LAGIS Trocar has applications in general, thoracic, gynecological, laparoscopic, abdominal and other minimally invasive surgical procedures to penetrate the body cavity and establish a path of entry for endoscopic instruments.
- Expected clinical benefits: reduce the surgery risks, shorten the surgery time, minimize the wound, and reduce the bleeding situation.
- Patient target group: the LAGIS Trocar is intended for adults who need minimally invasive procedures.

◆ **CONTRAINDICATIONS:**

The LAGIS Trocar is not intended for use when minimally invasive techniques are contraindicated.

◆ **INSTRUCTIONS FOR USE:**


Prior to use, inspect the contents of the packaging and the packaging itself to ensure the integrity is not compromised.

1. Insert the obturator into the cannula until the tip is exposed and make sure it is fully assembled.
2. Make an incision at the puncture site just long enough to accept the trocar to penetrate the body cavity.
3. Retract the obturator out of the cannula once the trocar is securely installed.
4. Perform the surgery according to the endoscopic surgical procedures.
5. Remove the cannula upon completion of the procedure.
6. Dispose of the used LAGIS Trocar properly and safely by healthcare professionals in accordance with the national waste management regulation.

◆ **WARNINGS AND PRECAUTIONS:**

- Do not use if the device or the sterile packaging is open or damaged, as this may cause an increase in the incidence of wound infections.
- The LAGIS Optical Trocar is designed to be used with rigid endoscopes. Do not use the LAGIS Optical Trocar with flexible endoscopes.
- If endoscopic instruments or accessories from different manufacturers are used together in a procedure, verify compatibility prior to performing the procedure.
- Using an endoscopic instrument with a diameter smaller than specified for the LAGIS Trocar may result in desufflation of the body cavity.

- The surgical site should be inspected for hemostasis after removing the LAGIS Trocar from the body cavity.
- This device is packaged and sterilized for single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocess, or resterilization may create contamination, infection, cross-infection, or even compromise the integrity of the device and lead to device failure which in turn may result in patient injury, illness or death.
- This device should be operated only by surgeons with adequate training and familiarity with surgical techniques. Consult medical literature for further information regarding techniques, hazards, contra-indications, and complications prior to performing the procedures.
- Return the device and packaging to the distributor or implement safe disposal handling by the healthcare professionals if the sterile packaging is damaged or unintentionally opened before use.
- The used device which comes into contact with bodily fluids should be considered as biomedical waste and safely disposed of by healthcare professionals in compliance with the national waste management regulation to prevent patients, users, and other persons from biological contamination or infection.
- In the event of malfunction of the device, changes in performance that they may affect safety, or even any serious incident, please report to the manufacturer, the competent authority of the Member State and the authorized representative immediately according to the contact information provided on the label or in the instructions for use.
- If the user feels the resistance is higher than usual when inserting the endoscopic instruments, the angle of insertion should be adjusted to align with the axis of the LAGIS Trocar for fear of damaging the valve.
(Please see the figure on the right-hand side.)
- Improper assembly or operation of the device may result in tissue or organ injuries from mechanical damage to the patient or operator.

◆ **ENVIRONMENTAL CONDITIONS FOR STORAGE:**

Appropriate storage environment is a clean and dry area away from sunlight with a temperature range of 13~30°C (55.4~86°F).

◆ **EXPLANATION OF SYMBOLS :**

	Temperature limit		Caution		Authorized representative in the European Community/European Union
	Do not re-sterilize		Manufacturer		Not made with natural rubber latex
	Do not re-use		Consult instructions for use or consult electronic instructions for use		Do not use if package is damaged and consult instructions for use
	Date of manufacture		Single sterile barrier system		Sterilized using ethylene oxide
	Medical device		Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.		

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