

www.lagis.com.tw/eifu_en.php

BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

◆ DEVICE DESCRIPTION :

The LAGIS Bipolar Coagulation and Cutting Forceps is a single use device sterilized by ethylene oxide. It is composed of forceps jaws, a cutting knife, a rotation knob, an insulating shaft, and a ratchet handle with a cutting trigger. The 5 mm diameter insulating shaft is designed for use through appropriate size trocars. With the ratchet handle, the forceps jaws are allowed to be locked in place. The device can be used in conjunction with compatible electrosurgical generators for electrocautery and coagulation.

PRODUCT MODEL

Code	Appropriate electrical generator
BFI-5330K-L	CONMED 2450
BFI-5330K-E	ERBE VIO 300D

◆ INDICATIONS/INTENDED PURPOSE :

- The LAGIS Bipolar Coagulation and Cutting Forceps is indicated for electrocautery, coagulation, mechanical cutting and grasping of tissue when connected to compatible electrosurgical generators during minimally invasive procedures.
- Expected clinical benefits: reduce the surgery risk, shorten the surgery time, minimize the wound, reduce the bleeding situation and easily controlled for the clinician.
- Patient target group: the LAGIS Bipolar Coagulation and Cutting Forceps is intended for adults who need endoscopic procedures.

◆ CONTRAINDICATIONS :

The LAGIS Bipolar Coagulation and Cutting Forceps is not intended for either contraceptive coagulation of fallopian tissue or for use when minimally invasive techniques are contraindicated; however, it may be used to achieve hemostasis following transection of the fallopian tube.

◆ INSTRUCTIONS FOR USE :

• BEFORE SURGERY

- Prior to use, inspect the contents of the packaging and the packaging itself to ensure the integrity is not compromised.
- Pre-test the device to verify the proper electrical functionality and generator setting:
 - Place the forceps jaws on a gauze pad that has been soaked with saline.
 - Connect the device to a compatible electrosurgical generator and make sure all connectors are securely connected and the electrosurgical generator's power setting is within the specified range (30-50W). The maximum voltage of the electrosurgical generators is 1500Vpk.
 - Activate the electrosurgical function by depressing the foot pedal.
 - Observe if steam is generated from the gauze pad, which indicates proper electrical performance.
- The LAGIS Bipolar Coagulation and Cutting Forceps can only be used when the abdominal cavity has been inflated and the laparoscope has been already set up. Insert the device through an appropriately sized trocar into the body cavity to perform endoscopic surgical operations.

• DURING SURGERY

- Ensure the electrosurgical generator power switch is OFF and the forceps jaws are closed prior to insertion and withdrawal through the trocar. Advance the LAGIS Bipolar Coagulation and Cutting Forceps through the trocar until it exits the trocar and is in view.
- To use the Ratchet mode, follow the steps below:
 - In the "Ratchet ON" mode: Ensure the forceps jaws are open and switch on the ratchet located on the handle, the forceps jaws will be locked and stay opened. Depress the lever for closure.
 - To release the forceps jaws, depress the lever again.
 - In the "Ratchet OFF" mode: The lever must be squeezed throughout the entire grasping, coagulation and transection process.
- To rotate the forceps jaws, rotate the rotation knob in either direction for the desired angle.
- Position the forceps jaws as desired for grasping, coagulation, and transection of tissue.
- Place the tissue between the forceps jaws and make sure it is placed in the center for the optimum electrosurgical function prior to depress the lever for closure. Do not grasp tissue for coagulation or transection beyond the serrated teeth. Ensure the forceps jaws are in contact only with the tissue to be grasped, coagulated, or transected.
- Activate the electrosurgical generator by depressing the foot pedal to coagulate until the coagulation effect has been achieved between the coagulating surfaces (approximately 3-6 seconds).

- To transect, depress the foot pedal to keep the electrosurgical generator activated. Advance the cutting knife slowly through the tissue by depressing the cutting trigger until the tissue has been transected.
- Discontinue the electrosurgical power application when coagulation or transection has been finished.
- After proper desiccation of tissue has been achieved, advance the cutting knife through the tissue by depressing the cutting trigger while maintaining a firm grasp of the tissue with forceps jaws.
- Upon completion of the procedures, retract the LAGIS Bipolar Coagulation and Cutting Forceps from the trocar with forceps jaws closed.

• AFTER SURGERY

- To clean the coagulating surfaces, use a gauze pad soaked with saline solution to clean the debris.
Note: Keep the cutting knife retracted while cleaning.
- Dispose the LAGIS Bipolar Coagulation and Cutting Forceps properly and safely by healthcare professionals in accordance with the national waste management regulation.

◆ WARNINGS AND PRECAUTIONS :

- Do not use if the device or the sterile packaging is open or damaged, as this may cause an increase in the incidence of wound infections.
- If trocars or accessories from different manufacturers are used together in the procedure, verify compatibility prior to performing the procedure and ensure that electrical insulation or grounding is not compromised.
- This device is packaged and sterilized for single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocess, or resterilization may create contamination, infection, cross-infection, or even compromise the integrity of the device and lead to device failure which in turn may result in patient injury, illness or death.
- This device should be operated only by surgeons with adequate training and familiarity with surgical techniques. Consult medical literature for further information regarding techniques, hazards, contra-indications, and complications prior to performing the procedures.
- Do not use power generators that may deliver more than 1500 Volts peak to the Endoscopic Instruments. The electrosurgical generator power level should be set as low as necessary to achieve the desired electrosurgical effect. Excessive power levels may result in instrument malfunction, insulation failure, and injury or thermal damage to patients or users.
- Always use the lowest possible output setting necessary to achieve the desired electrosurgical effect. In order to lessen the possibility of unintended burn injury, use the active electrode only for the minimum time necessary. Higher settings may result in spurious insulation failure alarms and/or insulation breakdown.
- Do not touch the coagulation surface while the device is activated. This may cause injury.
- The forceps jaws that are activated or hot from use can cause fire; do not place them near or in contact with flammable materials and substances (e.g. drapes, flammable gases, etc.)
- If forceps jaws come in contact with each other when the power is activated, instrument will short out. Release the forceps jaws so they are not in contact with each other and the device will become operational.
- Return the device and packaging to the distributor or implement safe disposal handling by the healthcare professionals if the sterile packaging is damaged or unintentionally opened before use.
- The used device which comes into contact with bodily fluids should be considered as biomedical waste and safely disposed of by healthcare professionals in compliance with the national waste management regulation to prevent patients, users, and other persons from biological contamination or infection.
- Do not attempt to repair the device. The use of a defective device or a device that becomes defective during surgery may result in loose parts that could drop into the field of operation.
- In the event of malfunction of the device, changes in performance that they may affect safety, or even any serious incident, please report to the manufacturer, the competent authority of the Member State and the authorized representative immediately according to the contact information provided on the label or in this instructions for use.
- Improper assembly or operation of the device may result in tissue or organ injuries from mechanical or thermal damage to the patient or operator.

◆ ENVIRONMENTAL CONDITIONS FOR STORAGE :

Appropriate storage environment is a clean and dry area away from sunlight with a temperature range of 13~30°C (55.4~86°F).

◆ EXPLANATION OF SYMBOLS :

	Temperature limit		Caution		Authorized representative in the European Community/European Union
	Do not resterilize		Manufacturer		Not made with natural rubber latex
	Do not re-use		Consult instructions for use or consult electronic instructions for use		Do not use if package is damaged and consult instructions for use
	Date of manufacture		Single sterile barrier system		Sterilized using ethylene oxide
	Medical device		Type BF applied part		Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	CE Marking according to MDR 2017/745				

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