
www.lagis.com.tw/eifu_en.php

BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

◆ DEVICE DESCRIPTION :

The LAGIS Wound Retractor is a single use device sterilized by gamma irradiation. It is composed of a green ring, a retraction sleeve, a blue ring, and an attachment ring. An incision template is included. Two elastic rings (the green ring and the blue ring) are located at both sides of the retraction sleeve, which is a one-piece design to perform incision retraction. In addition, the attachment ring helps to fix the retracted wound in a circular shape upon the wound retractor is installed.

◆ INDICATIONS/INTENDED PURPOSE :

- The LAGIS Wound Retractor is intended to provide maximum exposure with minimum incision size and to protect against wound contamination during the laparoscopic or open surgery.
- The expected clinical benefits: reduce surgical site infection and minimize incision length.
- Patient target group: The LAGIS Wound Retractor is intended for adults who need laparoscopic or open surgery.

◆ CONTRAINDICATIONS :

Use of the WOUND RETRACTOR is contraindicated when, in the judgment of the physician, use of the device would be contrary to the best interest of the patient.

PRODUCT MODEL

Models	Suggested Incisions Length	Cesarean Section	Compatible LAGIS Cap for wound retractor
WR-50-H	20-40 mm	NO	N/A
WR-60ES	20-40 mm	NO	CP-60ES
WR-60ES-H	20-40 mm	NO	CP-60ES
WR-80S	30-60 mm	NO	CP-80S
WR-110M	50-90 mm	NO	N/A
WR-110M-A	50-90 mm	NO	CP-110M
WR-160L	90-140 mm	YES	N/A
WR-190EL	110-170 mm	YES	N/A

◆ INSTRUCTIONS FOR USE :

Prior to use, inspect the contents of the packaging and the packaging itself to ensure the integrity is not compromised.

INSTALLATION AND ASSEMBLY :


1-Prepare the surgical site according to the hospital's procedure, making sure the skin is clean and dry.


2-Place the template over the incision site, and mark an appropriate incision line with a sterile marker (if needed).


3-Make an incision along the marked incision line.


4-Insert the blue ring through the incision.



5-Gently grasp the green ring at the opposite end and curl it inward.

6-Repeat step (5) until the wound is retracted and the wound retractor is securely installed.

7&8-Attach the attachment ring to the green ring to fix the retracted wound in a circular shape.

9-Perform the surgical procedure through the LAGIS Wound Retractor.

TO DISASSEMBLE :

1. Detach the attachment ring from the green ring.
2. Curl the green ring outward to loosen the wound retractor.
3. Remove the LAGIS Wound Retractor by grasping the blue ring out of the incision.

◆ **WARNINGS AND PRECAUTIONS :**

- Do not use if the device or the sterile packaging is open or damaged, as this may cause an increase in the incidence of wound infections.
- If trocars, endoscopic instruments or accessories from different manufacturers are used together in a procedure, verify compatibility prior to performing the procedure.
- This device is packaged and sterilized for single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocess, or resterilization may create contamination, infection, cross-infection, or even compromise the integrity of the device and lead to device failure which in turn may result in patient injury, illness or death.
- This device should be operated only by surgeons with adequate training and familiarity with surgical techniques. Consult medical literature for further information regarding techniques, hazards, contra-indications, and complications prior to performing the procedures.
- Return the device and packaging to the distributor or implement safe disposal handling by the healthcare professionals if the sterile packaging is damaged or unintentionally opened before use.
- The used device which comes into contact with bodily fluids should be considered as biomedical waste and safely disposed of by healthcare professionals in compliance with the national waste management regulation to prevent patients, users, and other persons from biological contamination or infection.
- In the event of malfunction of the device, changes in performance that they may affect safety, or even any serious incident, please report to the manufacturer, the competent authority of the Member State and the authorized representative immediately according to the contact information provided on the label or in this instructions for use.
- Improper assembly or operation of the device may result in tissue or organ injuries from mechanical damage to the patient or operator.

◆ **ENVIRONMENTAL CONDITIONS FOR STORAGE :**

Appropriate storage environment is a clean and dry area away from sunlight with a temperature range of 13~30°C (55.4~86°F).

◆ **EXPLANATION OF SYMBOLS :**

	Temperature limit		Caution		Authorized representative in the European Community/European Union
	Do not re-sterilize		Manufacturer		Not made with natural rubber latex
	Do not re-use		Consult instructions for use or consult electronic instructions for use		Do not use if package is damaged and consult instructions for use
	Date of manufacture		Single sterile barrier system		Sterilized using irradiation
	Medical device		Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.		
	CE Marking according to MDR 2017/745		UKCA Marking according to UK MDR 2002		

 **LAGIS[®] ENTERPRISE CO., LTD.**

No. 29, Gong 1st Rd., Dajia Dist., Taichung City 437, Taiwan
TEL: 886-4-26820767 FAX: 886-4-26820766
<http://www.lagis.com.tw>

 **Obelis s. a.**

Boulevard Général Wahis 53 1030 Brussels, BELGIUM
TEL: +(32)2.732.59.54 FAX: +(32)2.732.60.03
E-Mail: mail@obelis.net

Emergo Consulting (UK) Limited
c/o Cr360 - UL International Compass House, Vision
Park Histon Cambridge CB24 9BZ, United Kingdom

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