www.lagis.com.tw/eifu_en.php

BEFORE USING THE PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

◆ DEVICE DESCRIPTION:

The Tissue Containment Bag is a reinforced pouch with inner coating and drawstring closure.

◆ PRODUCT MODEL:

Product Code	Description	Total volume (ml)	Pouch diameter (mm)
EBS-0810	8*10 inch Specimen Bag	1500	130

◆ INDICATIONS/INTENDED PURPOSE:


- Tissue Containment Bag is applied to general, gynecological, thoracic, urological surgery, etc. for resection during the surgical procedure or tissues isolation before the surgery.
- Expected clinical benefits: avoid accidental spillage of content in the abdominal cavity.
- Patient target group: the Tissue Containment Bag is intended for adults who need laparoscopic surgical procedures.

◆ CONTRAINDICATIONS:


- Not intended for use when minimally invasive techniques are contraindicated.
- Not intended for use with any tissue that may not fit within the confines of the bag and allow complete closure of the bag.
- Not for use with power morcellators.

◆ INSTRUCTIONS FOR USE:

- First, fold over the corner of the Tissue Containment Bag, where a round hole is located, to form a small triangle.
- Use the grasping forceps to grasp the small triangle folded at the corner.
- Make sure that the grasping forceps have firmly grasped the Tissue Containment Bag, and then lock the handle of the grasping forceps.
- Roll the Tissue Containment Bag around the shaft for one full rotation in the direction of the folded triangle. Make sure the Tissue Containment Bag is firmly attached to the shaft of the grasping forceps first.
- Rotate the shaft along the direction of the folded triangle by turning the rotation knob on the handle of the grasping forceps.
- While turning the rotation knob, you can adjust the relative rotated position of the Tissue Containment Bag. Do NOT roll up the drawstring into the Tissue Containment Bag together.
- After the Tissue Containment Bag is rolled up into a tube shape, hold it tightly by hand.
- Then, insert the rolled Tissue Containment Bag into a trocar cannula of an appropriate size.
- Fix the rotation knob of the grasping forceps securely, to push the rolled Tissue Containment Bag completely into the trocar cannula through rotating the handle.
- Release the handle of the grasping forceps, and then take out the grasping forceps from the rolled Tissue Containment Bag.
- Use another grasping forceps as a support, to unfold the rolled Tissue Containment Bag.
- By using the two grasping forceps, place the specimen into Tissue Containment Bag.

After the specimen is placed into the Tissue Containment Bag, bind up the mouth of the bag by using the drawstring retained outside the trocar cannula.

Once the mouth of the Tissue Containment Bag is tightly closed by the drawstring, take out the trocar cannula first, and then retrieve the Tissue Containment Bag from the body cavity.

◆ WARNINGS AND PRECAUTIONS:

- The Tissue Containment Bag is not intended for use with sharp surgical instruments.
- Care should be taken to avoid damaging, puncturing or tearing the Tissue Containment Bag.
- Use caution when removing potentially malignant or harmful tissues, to isolate contamination from surrounding tissue.
- Do not use if the device or the sterile packaging is open or damaged, as this may cause an increase in the incidence of wound infections.
- This device is packaged and sterilized for single patient use only. Do not re-use, reprocess, or resterilize. Reuse, reprocessing, or resterilization may create contamination, infection, cross-infection, or even compromise the integrity of the device and lead to device failure which in turn may result in patient injury, illness or death.
- This device should be operated only by surgeons with adequate training and familiarity with surgical techniques. Consult medical literature for further information regarding techniques, hazards, contra-indications, and complication prior to performing the procedures.
- Return the device and packaging to the distributor or implement safe disposal handling by the healthcare professionals if the sterile packaging is damaged or unintentionally opened before use.
- The used device which comes into contact with bodily fluids should be considered as biomedical waste and safely disposed by healthcare professionals in comply with the national waste management regulation to prevent patients, users, and other persons from biological contamination or infection.
- In the event of malfunction of the device, changes in performance that they may affect safety, or even any serious incident, please report to the manufacturer, the competent authority of the Member State and the authorized representative immediately according to the contact information provided on the label or in this instructions for use.
- Improper assembly or operation of the device may result in prolong surgical time, patient allergy, and tissue or organ injuries from mechanical damage to the patient or operator.

◆ ENVIRONMENTAL CONDITIONS FOR TRANSPORT AND STORAGE:

Appropriate storage environment is a clean and dry area away from sunlight with a temperature range of 13~30°C (55.4~86°F).

◆ EXPLANATION OF SYMBOLS:

	Temperature limit		Caution		Authorized representative in the European Community/European Union
	Do not re-sterilize		Manufacturer		Not made with natural rubber latex
	Do not re-use		Consult instructions for use or consult electronic instructions for use		Do not use if package is damaged and consult instructions for use
	Date of manufacture		Single sterile barrier system		Sterilized using ethylene oxide
	Medical device		Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.		CE Marking according to MDR 2017/745

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