www.lagis.com.tw/eifu_en.php

BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

◆ DEVICE DESCRIPTION:

The LAGIS Insulated Suction Irrigation System consists of flexible tubing attached to a plastic housing/handle with a spring activated valve/stopper. Attached distally is a insulated stainless steel probe and proximally is a plastic housing cap. When properly connected and operated, the device irrigates and evacuates fluids during laparoscopic procedures. In the series of insulated stainless steel tubing, there is an insulating layer on the outside of the stainless steel tube. The function is to avoid transmitting electric current when used with other endoscopic electrode instruments.

Item	Pressure range	Flow rate
Value	0~1000 mmHg	>10L/min

PRODUCT MODEL

Models Numbers	Description	Probe Length(cm)
SA-450C-P	Direct Type Suction Irrigation Set for Dual Solution Bags	45
SA-360C-P		36
SA-450CT-P	Trumpet Type Suction Irrigation Set for Dual Solution Bags	45
SA-360CT-P		36
SI-0545P	Only Probe	45
SI-0536P		36

◆ INDICATIONS/INTENDED PURPOSE:

- The LAGIS Insulated Suction Irrigation System has applications in laparoscopic procedures to provide suction and irrigation functions to help flush blood and tissue debris from the operative site during laparoscopy to aid visualization.
- Expected clinical benefits: maintain clear view during laparoscopic procedures, assist in visualization of the bleeding point, remove blood pool, clear clots from the operative site and suction infected liquid.
- Patient target group: The LAGIS Insulated Suction Irrigation System is intended for adults who need laparoscopic procedures.


◆ CONTRAINDICATIONS:

The LAGIS Insulated Suction Irrigation System is not intended for use when endoscopic techniques are contraindicated.

◆ INSTRUCTIONS FOR USE:

Prior to use, inspect the contents of the packaging and the packaging itself to ensure the integrity is not compromised.

1. Pull the red guard tab located at the side of the housing(1) to relieve the buttons. (For Direct Type.)
2. For irrigation: insert the tubing spike(2) of The LAGIS Insulated Suction Irrigation System into a normal saline solution bag.
3. For suction: attach the suction tubing(3) of The LAGIS Insulated Suction Irrigation System to a suction canister.
4. Depress the blue button(4) to irrigate.
5. Depress the red button(5) to evacuate fluids.
6. Upon completion of the procedure, disconnect the tubes from normal saline solution bags and the fluid collecting canister.
7. Dispose of the used The LAGIS Insulated Suction Irrigation System, normal saline solution bags, and the canister properly and safely by healthcare professionals in accordance with the national waste management regulation.

◆ WARNINGS AND PRECAUTIONS:

- Do not use if the device or the sterile packaging is open or damaged, as this may cause an increase in the incidence of wound infections.
- This Instructions for use is to assist clinicians in the use of The LAGIS Insulated Suction Irrigation System. This is not an instructional reference for laparoscopic surgery.
- If trocars or accessories from different manufacturers are used together in a procedure, verify compatibility prior to performing the procedure.
- The probe and housing cap must be attached securely to prevent fluid leakage, as this may cause the device to lose function.
- This device is packaged and sterilized for single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocess, or resterilization may create contamination, infection, cross-infection, or even compromise the integrity of the device and lead to device failure which in turn may result in patient injury, illness or death.
- The normal saline solution only shall be used for irrigation solution.
- This device should be operated only by surgeons with adequate training and familiarity with surgical techniques. Consult medical literature for further information regarding techniques, hazards, contra-indications, and complications prior to performing the procedures.
- Return the device and packaging to the distributor or implement safe disposal handling by the healthcare professionals if the sterile packaging is damaged or unintentionally opened before use.
- The used device which comes into contact with bodily fluids should be considered as biomedical waste and safely dispose of by healthcare professionals in compliance with the national waste management regulation to prevent patients, users, and other persons from biological contamination or infection.
- In the event of malfunction of the device, changes in performance that they may affect safety, or even any serious incident, please report to the manufacturer, the competent authority of the Member State, and the authorized representative immediately according to the contact information provided on the label or in this instruction for use.
- Do not suck foreign objects or a large amount of tissue in a short period of time, as this may cause the device to lose suction function.
- Improper assembly or operation of the device may result in tissue or organ injuries from mechanical damage to the patient or operator.

Compatibility:

Please use LAGIS Insulated Suction Irrigation System with a normal saline solution bag and suction canister or irrigation/suction pump, tubing, 5 mm probe and 5 mm Trocar.

◆ ENVIRONMENTAL CONDITIONS FOR STORAGE:

Appropriate storage environment is a clean and dry area away from sunlight with a temperature range of 13~30°C (55.4~86°F).

	Temperature limit		Caution		Authorized representative in the European Community/European Union
	Do not re-sterilize		Manufacturer		Not made with natural rubber latex
	Do not re-use		Consult instructions for use or consult electronic instructions for use		Do not use if package is damaged and consult instructions for use
	Date of manufacture		Single sterile barrier system		Sterilized using ethylene oxide
	Medical device		Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.		CE Marking according to MDR 2017/745
	UKCA Marking according to UK MDR 2002				

ISI-0114 Ver. B 2025

LAGIS ENTERPRISE CO., LTD.

No. 29, Gong 1st Rd., Dajia Dist., Taichung City 437, Taiwan
TEL: 886-4-26820767 FAX: 886-4-26820766
http://www.lagis.com.tw

Obelis s. a.

Boulevard Général Wahis 53 1030 Brussels, BELGIUM
TEL: +(32)2. 732.59.54 FAX: +(32)2.732.60.03
E-Mail: mail@obelis.net

Emergo Consulting (UK) Limited
c/o Cr360 - UL International Compass House, Vision
Park Histon Cambridge CB24 9BZ, United Kingdom

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